

Rapydtest®

FOR THE DETECTION OF HUMAN HAEMOGLOBIN IN FAECES



FOB Rapydtest®



MICROBIOLOGY

SINGLE USE IN VITRO DIAGNOSTIC DEVICE

Intended Use

The Apacor FOB Rapydtest® is an immunochemical test device intended for the qualitative detection of faecal occult blood to be used in laboratories or physicians offices. It is a useful aid to detect bleeding caused by a number of gastrointestinal disorders, e.g., diverticulitis, colitis, polyps, and colorectal cancer. Faecal occult blood tests are recommended for use in:

1. routine physical examinations
2. routine hospital testing
3. screening for colorectal cancer or gastrointestinal bleeding from any source.

Explanation of the Test


Two types of FOB tests are commercially available: guaiac dye and immunochemistry. Guaiac tests are widely used but lack accuracy. Guaiac dye is a naturally occurring phenolic compound that can be oxidized to quinone by hydrogen peroxidase activity of human haemoglobin (hHb) resulting in a detectable colour change. The sensitivity and specificity of guaiac tests are much lower than those of immunochemical assays. The low accuracy of the guaiac dye tests is related to dietary peroxidases, including haemoglobin from meat and uncooked fruits and vegetables. Non-cancerous gastrointestinal tract bleeding and iron intake may also cause false positive results with guaiac tests¹.

Immunochemical tests are highly accurate for the detection of hHb compared to the guaiac method. The results of immunochemical FOB tests (iFOBT) are not affected by dietary peroxidases, animal blood or ascorbic acid. A Japanese study demonstrated that iFOBT screening tests reduced mortality of colorectal cancer by 60%². The Apacor FOB Rapydtest® is an iFOBT designed to specifically detect low levels of human faecal occult blood. It can be performed within 10 minutes by minimally skilled personnel and without the use of laboratory equipment.

Principle

The Apacor FOB Rapydtest® is a lateral flow chromatographic immunoassay. The test cassette consists of:

1. a burgundy coloured conjugate pad containing monoclonal anti-hHb antibodies conjugated with colloidal gold (anti-hHb conjugates) and a control antibody conjugated with colloidal gold
2. a nitrocellulose membrane strip containing a test band (T band) and a control band (C band). The T band is pre-coated with a monoclonal anti-hHb antibody, and the C band is pre-coated with a control band antibody.



A trace amount of haemoglobin (hHb) is first extracted from the faecal specimen with the sample extraction tube. When an adequate volume of test specimen is dispensed into the sample well of the test cassette, the specimen migrates by capillary action across the cassette. hHb, if present in the specimen at concentrations equal to or higher than 50ng/ml, will bind to the anti-hHb conjugates. The immunocomplex is then captured on the membrane by the pre-coated antibody forming a burgundy

coloured T band, indicating a positive test result. Absence of the T band suggests that the concentration of hHb in the specimen is below the detectable level, indicating a negative result.

The test contains an internal control (C band) which should exhibit a burgundy coloured band of the immunocomplex of control antibodies regardless of any colour development on the T band. If the C band does not develop, the test result is invalid and the specimen must be retested with another device.

Reagents and Materials Provided

1. Individually sealed foil pouches containing:
 - A. One cassette device
 - B. One desiccant
2. Sample extraction tubes, each containing 2ml extraction buffer
3. Patient ID stickers
4. One package insert (instruction for use)

Materials Required but not Provided

1. Clock or Timer
2. Faecal specimen container

Warnings and Precautions

For In Vitro Diagnostic Use

1. This package insert must be read completely before performing the test. Failure to follow the instructions may give inaccurate test results.
2. Do not open the sealed pouch, unless ready to conduct the assay.
3. Do not use any kit components beyond their stated expiration date.
4. Do not use the components from any other type of test kit as a substitute for the components from this kit.
5. Bring all reagents to room temperature (15°C-30°C) before use.
6. **Do not scoop faecal specimen as this may lead to excess faecal specimen that may block the sample well and result in an invalid test result.**
7. **Do not use specimen with visible blood for the testing.**
8. Wear protective clothing and disposable gloves while handling the kit reagents and clinical specimens. Wash hands thoroughly after performing the test.
9. Users of this test should follow 'Good Laboratory Practice' for biosafety.
10. Do not smoke, drink, or eat in areas where specimens or kit reagents are being handled.
11. Extraction buffer contains 0.3% Proclin-300. Avoid contact with skin or eyes. Do not ingest.
12. Dispose of all specimens and materials used to perform the test as biohazardous waste.
13. The test result should be read 5 to 10 minutes after a specimen is applied to the sample well of the device. Reading results after 10 minutes may give erroneous results.
14. Do not perform the test in a room with strong air flow, ie. electric fan or strong air-conditioning.

Reagent Preparation and Storage Instructions

Store unused or remaining test kits at 2-30°C. If stored at 2-8°C, ensure that the test kit is brought to room temperature (15-30°C) before opening. The test kit is stable through the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit to temperatures over 30°C.


Patient Preparation

1. A specimen should not be collected from a patient with following conditions that may interfere with the test results:
 - Menstrual bleeding
 - Bleeding haemorrhoids
 - Constipating bleeding
 - Urinary bleeding
2. Dietary restrictions are not necessary.
3. Alcohol and certain medications such as aspirin, indomethacin, phenylbutazone, reserpine, corticosteroids, and nonsteroidal anti-inflammatory drugs may cause gastrointestinal irritation and subsequent bleeding, and produce positive reactions. On the advice of a physician, these medicines may be temporarily discontinued for 7 days prior to and during the test period.

Specimen Collection and Handling

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.

1. Collect a stool specimen in a clean, dry receptacle.
2. Fill in all required information on the patient ID sticker and apply to the sample extraction tube.
3. Open the sample extraction tube by unscrewing the top and use the collection stick to randomly pierce the stool specimen in five different sites. Do not scoop stool specimen. Ensure that stool specimen is only in the grooves of the collection stick. Excess stool specimen may lead to an invalid test result.
4. Replace the collection stick and tighten securely to close the sample extraction tube.
5. Shake the stool collection device vigorously to extract the hHb in the specimen.




The specimen is now ready for testing, transportation or storage.

Note: It is recommended to test the specimen immediately after extraction. If not tested immediately, the extracted specimen may be stored at room temperature (20-37°C) for up to 10 days or at 2-8°C for up to 21 days. For longer storage, the extracted specimen may be frozen at -20°C. Avoid multiple freeze-thaw cycles.

Test Procedure

- STEP 1: Bring the specimen and test components to room temperature if refrigerated or frozen.
- STEP 2: When ready to test, open the pouch at the notch and remove the test cassette. Place the test cassette on a clean, flat surface.
- STEP 3: Shake the sample extraction tube vigorously to ensure an effective liquid suspension.
- STEP 4: Position the sample extraction tube upright and twist off the dispenser cap. Holding the sample extraction tube vertically, dispense 2 drops (70-90 µl) into the sample well of the cassette. Do not overload the sample well.



STEP 5: Start the timer.

STEP 6: Result can be read 5 -10 minutes after adding the sample. Positive results may be visible in a time period as short as 1 minute.

Do not read results after 10 minutes. To avoid confusion, discard the test device after interpreting the result.

Quality Control

1. **Internal Control:** This test contains a built-in control feature, the C band. The C band develops after adding specimen extract. If the C band does not develop, review the entire procedure and repeat the test with a new device.
2. **External Control:** Good Laboratory Practice recommends using external positive and negative controls to assure the proper performance of the assay, particularly under the following circumstances:
 - A - A New operator uses the kit, prior to performing testing of specimens.
 - B - A new lot of test kits is used.
 - C - A new shipment of kits is used.
 - D - The storage temperature of the kit falls outside of 2-30°C.
 - E - The temperature of the test area falls outside of 15-30°C.
 - F - To verify a higher than expected frequency of positive or negative results.
 - G - To investigate the cause of repeated invalid results.

Interpretation of Assay Result

1. **Negative Result:** If only the C band develops, the test indicates that the concentration of hHb in the specimen is below 50 ng/ml buffer. The result is negative or non-reactive.



2. **Positive Result:** If both the C and T bands develop, the test indicates that the concentration of hHb in the specimen is equal to or higher than 50 ng hHb/ml buffer. The result is positive or reactive.



3. **Invalid:** If no C band develops, the assay is invalid regardless of colour development on the T band as indicated below. Repeat the assay with a new device. **If caused by an excess amount of faecal specimen collected, collect a new specimen and retest.**



Performance Characteristics

SENSITIVITY

The analytical sensitivity of the test is 50 ng hHb/ml buffer or 7 mg hHb/g faeces approximately.

SPECIFICITY

The Apacor FOB Rapydtest® is specific for hHb. The following substances do not interfere with test results when spiked into negative specimens and weak positive specimens.

SUBSTANCE	CONCENTRATION	APACOR FOB ROPYDTEST®	
		NEGATIVE	POSITIVE
CONTROL	N/A	NEGATIVE	POSITIVE
BOVINE HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
CHICKEN HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
FISH HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
GOAT HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
HORSE HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
PORK HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
RABBIT HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
SHEEP HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
TURKEY HAEMOGLOBIN	2 mg/ml	NEGATIVE	POSITIVE
ASCORBIC ACID	20 mg/dl	NEGATIVE	POSITIVE
BILIRUBIN	100 mg/dl	NEGATIVE	POSITIVE
CAFFEINE	40 mg/dl	NEGATIVE	POSITIVE
GLUCOSE	2,000 mg/dl	NEGATIVE	POSITIVE
HORSERADISH PEROXIDASE	20 mg/ml	NEGATIVE	POSITIVE
IRON (Fe ³⁺ /Fe ²⁺)	5 mg/ml	NEGATIVE	POSITIVE

Aqueous extracts of cooked and uncooked beef, broccoli, cabbage, cantaloupe, cauliflower, chicken, fish, horseradish, lamb, parsnip, pork, red radish, turkey and turnip were spiked into negative and positive specimens. No interference was detected. Additionally, toilet bowl cleansers do not interfere with the results of the Apacor FOB Rapydtest®.

DOSE HOOK EFFECT

No false negative results due to the dose hook effect were observed for specimens containing hHb at concentrations up to 4 mg/ml.

Ordering Information

PRODUCT	PACK SIZE	CODE
FOB Rapydtest®	25	1642

Products can be ordered direct from Apacor or from an appointed distributor

Visit our website for all the latest information www.apacor.com or e-mail on: orders@apacor.com



UNIT 5, SAPPHIRE CENTRE,
FISHPONDS ROAD, WOKINGHAM,
BERKSHIRE, RG41 2QL, UNITED KINGDOM
TEL: +44 (0)118 979 5566
FAX: +44 (0)118 979 5186



MDSS GmbH
Schiffaraben 41
30175 Hanover
Germany

REPRODUCIBILITY

Known positive samples were tested in multiple assays and identically positive results were observed. Similarly, known negative samples produced negative results when tested in multiple assays.

CLINICAL PERFORMANCE

A total of 175 specimens were collected for the performance study. The Apacor FOB Rapydtest® was compared to a leading commercial iFOB Rapid Test marketed in the U.S. and other regions. Comparison for all specimens is shown in the following table:

REFERENCE TEST	APACOR FOB ROPYDTEST®		
	POSITIVE	NEGATIVE	TOTAL
POSITIVE	47	1	48
NEGATIVE	1	126	127
TOTAL	48	127	175

Relative Sensitivity: 97.9%, Relative Specificity: 99.2%, Overall Agreement: 98.9%

Limitations of the Test

1. Test Procedure and the Interpretation of Assay Results must be followed closely when testing for the presence of occult blood in faeces. Failure to follow the procedure may give inaccurate results.
2. The Apacor FOB Rapydtest® is to aid in diagnosis and is not intended to replace other diagnostic procedures such as G.I. fibroscope, endoscopy, colonoscopy, or X-ray analysis. Test results should not be deemed conclusive with respect to the presence or absence of gastrointestinal bleeding or pathology. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source for the occult blood in the faeces.
3. A negative or non-reactive result can be obtained even when a gastrointestinal disorder is present. For example, some polyps and colorectal cancers may bleed intermittently or not at all during certain stages of the disease. A negative or non-reactive result can also be obtained if the quantity of occult blood present in the specimen is below the detection limit of the assay.
4. The Apacor FOB Rapydtest® has not been validated for testing of patients with hemaglobinopathies.
5. Specimens containing visible blood may produce negative results due to the hook effect.

References

1. Allison JB, Takawa IS, Ransom LJ, Adrian AL. A comparison of fecal occult blood tests for colorectal -cancer screening. N. Eng. J. Med. 1996; 334:155-159.
2. Saito H. Screening for colorectal cancer by immunochemical fecal occult blood testing (Review). Jpn J. Cancer Res 1996; 87:1011-1024.

SAFETY DATASHEET

FOB ROPYDTEST®



This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

SECTION 1 IDENTIFICATION OF THE SUBSTANCE/ MIXTURE AND THE COMPANY/UNDERTAKING

1.1 Product Identifier: 1642 FOB Rapydtest® Sample Extraction Buffer

1.2 Relevant identified uses of the substance or mixture and uses advised against: Sample Extraction Buffer for use with the above Rapydtest® in hospitals and laboratories. For *in vitro* use only.

1.3 Details of the supplier of the Safety Data Sheet:
Apacor Limited, Unit 5 Sapphire Centre, Fishponds Road, Wokingham, Berkshire, RG41 2QL, United Kingdom
+44 (0) 118 979 5566
technical@apacor.com

1.4 Emergency telephone number:
+44 (0)118 979 5566
(Monday-Friday 0900-1700 excluding UK Public Holidays)

SECTION 2 HAZARDS IDENTIFICATION

2.1 Classification of the substance or mixture
Classification according to Regulation (EC) No 1272/2008 [CLP]:
Not classified as hazardous in concentration of <0.6% (Proclin-300).

2.2 Label elements
Labelling according to Regulation (EC) No 1272/2008 [CLP]

The product does not contain a hazardous ingredient in an amount that requires identification and labelling according to the concentration limit/cut-off values of EC directives. This product contains no hazardous constituents, or the concentration of all chemical constituents is below the regulatory threshold limits described by Occupational Safety Health Administration Hazard Communication Standard 29 CFR 1910.1200 and the European Directive 91/155/EEC, 93/112/EC and (EC) 1272/2008 (CLP).

Pictogram: None

Signal Word: -

Hazard statement(s): -

Precautionary statements:-

2.3 Other hazards

Bio-hazards: All the biological substances are derived from *in vitro* culture system or animal materials which are free of known-pathogens for human. Thus, no bio-hazards can be claimed in the product.

SECTION 3 COMPOSITION/INFORMATION ON INGREDIENTS

3.2 Mixtures

Hazardous ingredients according to Regulation (EC) No 1272/2008

Component: **Proclin 300**

≤ 0.5% Proclin 300 [0.0015% or 0.015% active ingredients – reaction mass of: 5- chloro-2-methyl-4-isothiazolin-3-one (C₄H₄ClNOS; CAS# 26172-55-4, EC No 247-500-7) and 2-methyl-2H -isothiazol-3-one (C₄H₅NOS; CAS# 2682-20-4, EC No 220-239-6) (3:1)].

CAS No: 55965-84-9

EC No: 613-167-00-5

Classification: Skin Sens. Cat 1 (H317); Skin Irrit. Cat 2 (H315); Eye Irrit. Cat 2 (H319)

Concentration: < 0.3%

SECTION 4 FIRST AID MEASURES

4.1 Description of first aid measures

Consult a physician. Show this safety data sheet to the doctor in attendance.

If inhaled: Inhalation of any component in this kit is unlikely. If a component of this kit is inhaled and causes discomfort, move exposed individual to fresh air. Seek medical attention if breathing is difficult or symptoms persist.

In case of skin contact: The Sample Extraction Buffer is hazardous by skin contact. In case of contact, immediately clean skin with plenty of water. Remove contaminated clothing and shoes. Cold water may be used. Wash clothing before reuse. Thoroughly clean shoes before reuse. The animal proteins and dried reagents absorbed into the nitrocellulose membrane and the fibre conjugate pad are very unlikely to be hazardous by skin contact, but cleaning the skin after use is advisable.

If swallowed: Ingestion of small amounts of the Sample Extraction Buffer is toxic; a physician should be consulted immediately. The animal proteins and dried reagents absorbed into the nitrocellulose membrane and the fibre conjugate pad are very unlikely to be ingested or be hazardous by ingestion. However, a physician should be consulted should ingestion occur.

In case of eye contact: The test device is very unlikely to come into contact with the eye, however, a physician should be consulted should contact occur. In case of contact with the Sample Extraction Buffer, immediately flush eyes with plenty of water for at least 15 minutes. Cold water may be used. Get medical attention.

Aggravating Condition: Repeated or prolonged exposure is not known to aggravate medical conditions.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (Section 2.2) and/or Section 11.

4.3 Indication of any immediate medical attention and special treatment needed

No data available.

SECTION 5 FIRE FIGHTING MEASURES

5.1 Extinguishing media

For small fires, use dry chemical, carbon dioxide, or alcohol-resistant foam. No direct contact with water.

5.2 Special hazards arising from the substance or mixture

When involved in a fire, this material can decompose and produce irritating fumes and toxic gases (e.g., carbon monoxide, carbon dioxide, sulfuric dioxide).

SAFETY DATASHEET

FOB ROPYDTEST®



This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

Explosion Sensitivity to Mechanical Impact: Not sensitive under normal conditions.

Explosion Sensitivity to Static Discharge: Not sensitive under normal conditions.

5.3 Advice for firefighters

This material will not significantly contribute to the intensity of a fire. Use extinguishing material suitable to the surrounding fire. Utilize proper personal protective equipment when responding to any fire. Incipient fire responders should wear eye protection. Structural firefighters must wear Self-Contained Breathing Apparatus and full protective equipment. Move containers from fire area if it can be done without risk to personnel. If possible, prevent runoff water from entering storm drains, bodies of water, or other environmentally sensitive areas.

SECTION 6 ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures

Lab coat and gloves.

6.2 Environmental precautions

No data available.

6.3 Methods and material for containment and cleaning up

Use absorbent paper towel or cloth to absorb the spill solution and dispose or clean the contaminated surface in accordance with local procedures or appropriate standards

6.4 Reference to other sections

For disposal, see Section 13.

SECTION 7 HANDLING AND STORAGE

7.1 Precautions for safe handling

Do not eat, drink, smoke or apply cosmetics in laboratory area. Use the product according to the product insert.

7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed. Keep product at 2-30°C. Do not freeze or expose to temperature higher than 30°C. Keep away from children.

7.3 Specific end use(s)

No other specific uses are specified apart from those listed in Section 1.2.

SECTION 8 EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

CAS #	Chemical Name	OSHA (PEL)	ACGIH (TLV)*	NIOSH
26172-55-4	Proclin 300	N/A	0.1mg/m ³	N/A

Biological Exposure Index (ACGIH).

Other exposure limits for potential decomposition products:
None.

8.2 Exposure controls

Engineering Control: Eye bath. Use adequate ventilation to keep airborne concentrations low.

Hygiene Measures: Wash hands after handling compounds and before eating, smoking, using lavatory, and at the end of the day.

Personal Protective Equipment:

Respiratory Protection: None needed under normal conditions of use.

Skin and Body: Lab coat as indicated by general lab practice guidelines.

Eyes: Safety glasses or face shield are recommended to prevent eye contact.

Hand: Compatible chemical resistant gloves.

SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

a) **Appearance** Form: liquid

b) **Odour** odourless

c) **Odour threshold** no data available

d) **pH** 4.1

e) **Melting point / freezing point** no data available

f) **Initial boiling point and boiling range** 233°C

g) **Flash point** non-combustible

h) **Evaporation rate** <1

i) **Flammability (solid, gas)** no data available

j) **Upper/lower flammability or explosive limits** not applicable

k) **Vapour pressure** 0.06mmHg

l) **Vapour density** not applicable

m) **Relative density** no data available

n) **Solubility (ies)** soluble

o) **Partition coefficient: n-octanol/water** no data available

p) **Auto-ignition temperature** not applicable

q) **Decomposition temperature** no data available

r) **Viscosity** no data available

s) **Explosive properties** no data available

t) **Oxidising properties** no data available

9.2 Other information

No data available

SECTION 10 STABILITY AND REACTIVITY

10.1 Reactivity

No data available.

10.2 Chemical stability

Stable under normal storage conditions.

10.3 Possibility of hazardous reactions

Strong oxidizers, Perchloric Acid, > 25°C for prolonged periods.

10.4 Conditions to avoid

Contact with strong oxidizers, Perchloric Acid, > 25°C for prolonged periods.

SAFETY DATASHEET

FOB ROPYDTEST®



This Safety Datasheet complies with the requirements of Regulation (EC) No 1907/2006

10.5 Incompatible materials

Strong oxidizers, Perchloric Acid, > 25°C for prolonged periods.

10.6 Hazardous decomposition products

Hydrogen chloride, oxides of nitrogen and sulphur.

SECTION 11 TOXICOLOGICAL INFORMATION

11.1 Information of toxicological effects

No adverse effects on the health are expected from the components of the product. There is no aquatic toxicity data for this product at this time. Individual aquatic toxicity studies have been completed for the below listed chemicals.

Proclin 300

Acute toxicity: no data available

Skin corrosion/irritation: no data available

Serious eye damage/eye irritation: no data available

Respiratory or skin sensitisation: no data available

Germ cell mutagenicity: no data available

Carcinogenicity: no data available

Reproductive toxicity: no data available

Specific target organ toxicity - single exposure: no data available

Specific target organ toxicity - repeated exposure: no data available

Aspiration hazard: no data available

Additional Information

Proclin 300	LD50 oral rat, female: 3723 mg/kg
	LD50 oral rat, male: 3600 mg/kg
	LD50 skin rabbit, female: >3600 mg/kg
	LD50 skin rabbit, male: 3500 mg/kg

To the best of our knowledge, the chemical, physical and toxicological properties of this chemical have not been thoroughly investigated.

SECTION 12 ECOLOGICAL INFORMATION

12.1 Toxicity

No adverse effects on the environment are expected from the components of this kit. There is no aquatic toxicity data for this kit at this time.

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

There is limited potential for the components within this kit to accumulate in plant or animal systems.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted.

12.6 Other adverse effects

No data available.

12.7 Additional information

None.

SECTION 13 DISPOSAL CONSIDERATIONS

13.1 Waste treatment methods

Waste must be disposed of in accordance with federal, state and local environmental control regulations.

SECTION 14 TRANSPORT INFORMATION

14.1 UN number none

14.2 UN proper shipping name none

14.3 Transport hazard class(es) This substance is considered to be non-hazardous for transport.

14.4 Packing group none

14.5 Environmental hazards Do not discharge effluent containing this kit into streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge.

Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance, contact the appropriate environmental agency.

14.6 Special precautions for user no data available

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Not intended to be transported in bulk.

SECTION 15 REGULATORY INFORMATION

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

No data available.

15.2 Chemical Safety Assessment

No chemical safety assessment has been carried out for this product.

SECTION 16 OTHER INFORMATION

Full text of H-Statements referred to in Section 3

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

Eye Irrit. Eye irritation.

Skin Sens. Skin sensitisation.

Skin Irrit. Skin irritation.

The information supplied in this SDS is correct to the best of our knowledge. We do not accept any liability for loss, injury or damage, which may result from its use.



MDSS GmbH
Schiffaraben 41
30175 Hanover
Germany