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## SOLVENT FREE FAECAL PARASITE CONCENTRATOR

APACOR


### Health and Safety Benefits

- Totally enclosed/sealed process
- Solvent free protocol
- Disposable device
- Single use, no sample contamination
- Ready to use systems available

### Performance Benefits

- Optimum sample recovery
- Enhanced sample clarity
- Rapid four step process
- Human resources optimised
- Reduce human error

Code	Pack Size
108000	1000
108100	1000
108800	50

### Patented Filter

Large particles are rejected without obscuring the filtration.

Recovery rate with Parasep<sup>®</sup> is comparable to traditional sieve method, ie: Ridley-Allen.

### Debris Trap

Rejected particles are trapped to prevent extrusion into the Sedimentation Cone during centrifugation.

### Air/Liquid Seal and Safety Lock

The 'seal' prevents the release of biohazardous material. The 'lock' ensures the Mixing Chamber and Filter are removed together for safe disposal.

### Sedimentation Cone

Parasites are efficiently collected in the sedimentation cone during centrifugation and easily aspirated/pipetted for microscopic examination.

### Centrifuge Compatability

Designed to fit all 15ml centrifuge buckets.

### Procedure

See label for storage conditions and expiry date. Please adhere to the following guidelines when handling Mini Parasep<sup>®</sup> SF. To avoid cross contamination, the Mini Parasep<sup>®</sup> SF device should remain closed at all times except when introducing the sample or when retrieving the final concentrated sample for examination.

#### STEP 1 - SAMPLE PREPARATION

**Preserved Samples**  
Shake or vortex the incoming preserved sample to thoroughly mix.  
Transfer either 2ml (**US Gold Standard**) or 3ml (**ARUP JCM**) of the emulsified stool into the Mini Parasep<sup>®</sup> SF mixing chamber. In the event of:  
**Thick Stool Samples**—please add 10 drops of Apacor Triton X, then please enclose and vortex/shake to emulsify prior to transferring the sample;  
**Liquid Stool Samples**—please add 4ml instead of 2ml or 3ml to ensure that a sediment is formed after centrifugation is performed.

#### STEP 2 - EMULSIFICATION

Seal the Mini Parasep<sup>®</sup> SF by screwing in the filter/sedimentation cone unit.

#### STEP 3 - CENTRIFUGATION

Invert Mini Parasep<sup>®</sup> SF and perform centrifugation at:  
**500g** for ten minutes (**US Gold Standard**);  
**400g** for two minutes (**ARUP JCM**).

NOTE: TO CALCULATE THE REQUIRED RPM FOR ANY CENTRIFUGE.

$$RPM = \sqrt{\frac{g}{1.12r}} \times 1000$$

RPM - rotor speed in revs/min.  
g - centrifugal force (max.1000g)  
r - radius, horizontal distance between sedimentation cone tip and spindle centre measured in mm.


#### STEP 4 - EXAMINATION

Unscrew and discard the filter and mixing tube.

Decant the supernatant.

Transfer sediment to slide to perform examination where all temporary (Lugol's Iodine) and permanent staining techniques (Trichrome) can be conducted providing a universal fixative is used.

If a universal fixative is used in conjunction with this device, molecular and EIA tests can be conducted from the sediment.

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