
# Leishmania Rapydtest®

## Instructions for Use

### EN

#### Specimen Collection

The Apacor Leishmania Rapydtest® can be performed on either serum or plasma. Remove the serum or plasma from the clot of red cells as soon as possible to avoid haemolysis. Only clear, non-haemolysed specimens should be used.

Testing should be performed as soon as possible after sample collection. Do not leave samples at room temperature for prolonged periods. Specimens can be refrigerated at 2-8°C for up to 3 days. For longer periods specimens should be stored below -20°C. Bring specimens to room temperature prior to testing. The frozen specimens must be completely thawed prior to testing. Specimens should not be repeatedly frozen and thawed. If specimens are to be shipped, they should be packed in compliance with regulations covering the transportation of etiologic agents.

#### Test Procedure

Read the entire procedure carefully prior to performing any tests. Allow serum samples to reach room temperature (20-30°C) prior to testing. Remove the dipstick from the foil pouch. Use the test as soon as possible but within 1 hour after removal from pouch to prevent exposure to humidity.

Add 20µl of serum or plasma to the dipstick in the area below the arrows.

Place the dipstick into a test tube, or well of a 96 well tissue culture plate so that the end of the strip is facing downward as indicated by the arrows on the strip.

Add 2-3 drops (150µl) of the Chase Buffer solution provided with this test kit to the base of the well or test tube.

Read the results in 10 minutes. It is significant that the background is clear before reading the test, especially when sera have low titer of anti-leishmanial antibody and only a weak band appears in the test region (T). Results interpreted after 10 minutes can be misleading.

#### Quality Control

A procedural control is included in the test. A red line appearing in the control region is considered as an internal positive procedural control. A clear background in the control region is considered as an internal negative procedural control.

#### Interpretation of Results

**POSITIVE:** Two distinct red lines will appear, one in the test region and another in the control region.

**NEGATIVE:** Only a single red line appears in the control region. No apparent red or pink line appears in the test region.

**INVALID:** No lines appear at either the control or test line area. The test is also invalid if no control line appears, but a test line is seen. It is recommended to retest using a new dipstick and fresh serum. If the problem persists, discontinue using the lot immediately and contact Apacor.

**NOTES:** The red colour in the test line region will vary depending on the concentration of anti-leishmanial antibodies present. However, neither the quantitative value nor the rate of increase in antibodies can be determined by this qualitative test.

#### Limitation of the Procedure

The Apacor Leishmania Rapydtest® is for in vitro diagnostic use only. This test should be used for the detection of antibodies to Leishmania in serum or plasma samples. This test will only indicate the presence of antibodies to Leishmania the specimen and should not be used as the sole criterion for the diagnosis of Leishmania. As with all diagnostic tests, all results must be considered with other clinical information available to the physician. Persons with advanced HIV infection or other immunocompromised diseases frequently have low or undetectable anti-Leishmania antibodies. If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is recommended. A negative result at any time does not preclude the possibility of Leishmania.

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Apacor Limited
MDSS GmbH