ParaClick®

REF 250000

Instructions for use

INTENDED USE

The ParaClick® pinworm paddle is a sampling and microscopy device for the detection of enterobiasis ova to aid in the diagnosis of enterobiasis infection.

SUMMARY

The human pinworm, scientifically known as *Enterobius vermicularis* (*E. vermicularis*), is a nematode (round-worm) which is estimated to infect 200 million people worldwide with children aged 5~10 years old accounting for over 30% of cases (Burkhart and Burkhart., 2005).

E. vermicularis, commonly known as pinworm due to the female's long pointed tail, reside in the cecum, appendix, and ascending colon. Female pinworms are 8-13 mm long, and males are 2-5 mm long. Enterobiasis may be asymptomatic or may present clinical manifestations in the form of perianal pruritus due to female worms depositing eggs in the perianal region leading to local epidermis irritation and secondary bacterial infections, abdominal discomfort, loss of appetite, weight loss, insomnia, restlessness, irritability (Vermund & Wilson, 2000; Burkhart & Burkhart, 2005). The highly contagious nature of this disease may be due to faecal-oral route of transmission. In addition, contamination of bed linens, clothes and hands may play a significant role in transmission. After ingestion of eggs, it takes around 4 to 8 weeks for the appearance of new eggs around the anus.

The diagnosis of pinworm eggs and worms can be tricky in some individuals due to low parasite burden and no symptoms. Examining stool samples for pinworms is not recommended as these worms are often sparse in stool. Moreover, serologic tests are not available for diagnosing pinworm infections. Because of low invasiveness of the worms, neither blood eosinophilia nor elevated immunoglobulin E levels are expected to change. Therefore, it is recommended for microscopic identification of ova to be used for successful diagnosis (Wendt et al., 2019)

The ParaClick Pinworm Paddle is an all-in-one sampling and microscopy device for the detection of enterobiasis ova. The test uses an adhesive patch to collect ova, before folding to place the ova under the viewing

window, ready for microscopic analysis.

PRINCIPLE

The ParaClick® pinworm paddle is a collection and microscopy device for the detection of pinworm ova in suspected pinworm infection.

In this test, the adhesive section can be exposed, before being pressed against the perianal area to collect any potential pinworm ova. Once the sample has been collected, the adhesive area can be folded to sit under the built-in viewing window, allowing for microscopic examination with no further preparation. Presence of ova when viewed under a microscope is indicative of a positive result.

REAGENTS

Not applicable

MATERIALS

Provided: ParaClick Pinworm Paddles, Patient Insert (visual aid), Plastic Bags.

Not Provided: Microscope.

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- Do not use test if damaged.
- Do not eat, drink, or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves, and
- eye protection when specimens are being evaluated.
- The used test should be discarded according to local regulations.
- Humidity and temperature can adversely affect results.

STORAGE AND STABILITY

- Store as packaged at room temperature.
- Store the used tests in the refrigerator.
- The test is stable through the expiration date printed on the box.
- Do not use beyond the expiration date.

SPECIMENS

1. Pinworm ova detection is improved by collecting a sample on three consecutive days, as repeated samples are more likely to lead to successful collection of ova.
2. Ova detection is also more likely if samples are collected first thing in the morning, prior to a bowel movement and prior to washing.
3. Collected specimens should be stored in the fridge between collection and microscopy.

DIRECTIONS FOR USE / Procedure

Perform the test first thing in the morning, prior to bowel movements or washing.

Steps for patient:

1. Ensure the child's anal area is clean before going to bed.
2. Do the test first thing in the morning immediately after the child wakes (before a bowel movement and before washing).
3. Label one of the sample bags clearly with the date of sample collection using a ballpoint pen or permanent marker.
4. Wash your hands with soap and warm water and dry well on a clean towel or paper towel.
5. Remove one of the tests from the cardboard holder.
6. Peel the protective film from the adhesive area- DO NOT TOUCH THE ADHESIVE AREA.
7. Gently press the adhesive area several times on the skin around the perianal area.
8. Fold the test to seal and protect the sample. The test will 'click' shut.
9. Place the used test in the plastic bag previously

labelled with the collection date, seal the bag and place in the refrigerator.

10. Wash your hands.
11. Store the used tests in the refrigerator until all three samples have been collected.
12. After taking the third test, place all three samples in the envelope provided and return to your health care provider.

*Refer to ParaClick patient visual aid APA-VA-01 (English & Spanish)

Steps for microscopist:

1. Place the closed ParaClick onto the microscope stage, with the lens over the dipped viewing window.
2. Adjust the focus and light source for optimum/ comfortable viewing.
3. If present, enterobiasis ova should be visible at 10x magnification.
4. Once examined, dispose of the ParaClick as per local regulations

INTERPRETATION OF RESULTS

1. No pinworm eggs (*Enterobius vermicularis*) identified.
2. Positives reported as few, moderate, or many eggs identified.

QUALITY CONTROL

1. Implement internal quality control measures within the laboratory setting to monitor the accuracy and precision of the testing process.
2. Regularly check and calibrate laboratory equipment, such as microscopes.


LIMITATIONS

1. The ParaClick Pinworm Paddle is for in vitro diagnostic use only. The test should be used for the detection of human enterobiasis ova only.
2. As with all diagnostic tests, all results must be interpreted together
3. with other clinical information available to the physician.
4. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not

at any time preclude the possibility of enterobiasis infection with low concentration of ova.



 Apacor Limited
Sapphire Centre
Fishponds Road
Wokingham
Berkshire
RG41 2QL
United Kingdom

www.apacor.com
+44 (0) 1118 979 5566

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www.mdss.com
+49 511 6262 8630

	= CE marking (European Regulation 2017/746-CE on in vitro diagnostic medical devices)
	= Use by
	= Lot number
	= Catalogue number
	= Manufacturer
	= Contains sufficient for <n> tests
	= <i>In vitro</i> diagnostic medical device
	= Temperature limitation
	= Consult instructions for use
	= European Authorised Representative
	= Date of manufacture
	= Do not use if package is damaged
	= Caution

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MDSS GmbH