


RAPYDTEST® USING

COMBINED HRP2 / PLDH TECHNOLOGY FOR THE DIFFERENTIAL
DIAGNOSIS OF PLASMODIUM FALCIPARUM AND THE OTHER
PLASMODIUM SPECIES



APACOR

Rapydtest®
CareUS™
Malaria



PARASITOLOGY

SINGLE USE IN VITRO DIAGNOSTIC DEVICE

Performance Benefits

- Isolates Plasmodium falciparum HRP2 and Pan
- Specific LDH (Pf, Pv, Po, Pm) on separate test lines
- Combined antigen technology gives you increased accuracy
- User friendly cartridge format for ease of use and storage
- Integral vents prevent sample 'back flow' interference
- Results in 20 minutes

Intended Use

For the rapid qualitative detection of malaria HRP2 (histidine-rich protein 2) of *P. falciparum* and pLDH (plasmodium lactate dehydrogenase) of *P.falciparum*, *P.vivax*, *P.ovale* and *P.malariae* in human blood as an aid in the diagnosis of malaria infection.

Summary

Malaria is a serious parasitic disease characterized by fever, chills, and anemia. It is caused by a bite of infected Anopheles mosquitoes resulting transmission of protozoan parasite to human. There are four types of human malaria: *Plasmodium falciparum*, *P. vivax*, *P. ovale*, and *P. malariae*. When infected, the parasites (called sporozoites) migrated to the liver where they mature and then they are released to the blood stream of human infecting the red blood cells. Malaria infection occurs in more than 90 countries worldwide, but is mostly prevalent in sub-Saharan Africa. It is estimated that there are over 500 million clinical cases and 2.7 million malaria-caused deaths per year.

The careUS™ Malaria Combo Pf/PAN (HRP2/pLDH) Ag contains a membrane strip, which is pre-coated with two monoclonal antibodies as two separate lines across the test strip. One monoclonal antibody (test line 2) is PAN specific to pLDH of the Plasmodium species (*P.falciparum*, *P.vivax*, *P.malariae*, and *P.ovale*) and the other line (test line 1) consists of a monoclonal antibody specific to HRP2 of the *P.falciparum*. The conjugate pad is dispensed with antibodies absorbed on gold particles, which are specific to pLDH of PAN and specific to HRP2 of *P. falciparum*.

The careUS™ Malaria Combo Pf/PAN (HRP2/pLDH) Ag is designed for the differentiated diagnosis of *P.falciparum*, *P.vivax*, *P.malariae*, and *P.ovale* infection.

Precautions

In order to obtain reproducible results, the following rules must be observed:

- Test kits are for in vitro diagnostic and professional use only.
- Read the provided instructions for use before using the test kit and follow the provided information when using the test kit.
- Wash the hands thoroughly before and after using the test kit.
- Wear the protective gloves at all times while using the test kits and dispose the gloves immediately after each testing.
- Change the gloves and wash hands when contacted with potentially infectious materials.
- Do routine clean-up using an appropriate disinfectant.
- Do not eat or smoke while using the test.

- Use the test cassette and accessory components (lancet and alcohol swab) immediately after opening its package.
- All provided materials are single-use. Do not re-use any of the materials.
- Do not use the materials from different lots.
- Do not use a test cassette, if the packaging of the test cassette is compromised.
- Do not use the test kit if the expiration date has past.
- Lancet is sterile. If a lancet cap is loose and/or damaged, do not use the lancet.
- Do not swallow the assay buffer solution.
- Do not use the alcohol swab if the package is damaged.
- Observe the storage condition indicated on the packaging and box.
- Do not freeze the test kit. The refrigerated test cassettes need to be equilibrated to the room temperature prior to use.
- Dispose of waste in accordance with the local regulations.
- The buffer bottle should be lightly closed after each use and stored in a cool area avoiding direct sun light. Discard the buffer solution and bottle after the expiration date indicated on the bottle.
- The test may produce a false positive result for a patient with acute schistosomiasis.

Storage & Stability

The sealed pouch containing the test strip is designed to be stored at 1°C - 40°C for the duration of its shelf life. The bottle containing the Assay Buffer is designed to be stored at 1°C - 40°C for the duration of its shelf life. Exposure to temperatures over 40°C can impact the performance of the test and should be minimized. The strips should not be frozen. The test should be used within 15 minutes after removal from the pouch to prevent exposure to humidity. This product is stable for 30 months.

Specimen collection and storage

1. Clean area to be lanced with an alcohol swab.
2. Squeeze the end of the fingertip and pierce with a sterile lancet provided.
3. Wipe away the first drop of blood with sterile gauze or cotton.
4. Take a sample pipette provided, and collect the blood sample (5µl).

Accessories supplied

Sample pipette


Disposable capillary dropper



Capped Lancet






Alcohol pad






Test Procedure


1. Add 5µl of whole blood into the sample well (small well).
2. Add **two drops (60µl)** of Assay Buffer into the buffer well (updated procedure).
3. Read test result in 20 minutes




NEGATIVE




INVALID



POSITIVE
P.falciparum



POSITIVE
P.vivax, P.malariae
or *P.ovale*



POSITIVE
P.falciparum
or mixed infection

Results

Interpretation of the test

1. **Negative reaction**
The presence of only one band in the Control Area within the result window indicates a negative result.
2. **Invalid**
The test is invalid if the line in the Control Area does not appear. If this occurs, the test should be repeated using a new strip.
3. **Positive reaction - *P. falciparum***
The presence of three colour bands (three bands in the Control, "2" and "1" areas) or two bands (one band in the Control Area and another band in the "1" area) indicates a positive result for *P. falciparum*.
4. **Positive reaction**
P. vivax, P. malariae, or P. ovale. The presence of two colour bands (one band in the Control Area and another band in the "2" area) indicates a positive result for *P. vivax, P. malariae, or P. ovale.* The pLDH present in the sample reacts with the pan anti-pLDH conjugate and move through the test strip where the pLDH is captured by pan specific anti-pLDH.
5. **Positive reaction - mixed infection of *P. falciparum* and other species**
The presence of three colour bands (bands in the Control, "2" and "1" areas) indicates a positive result for *falciparum* or mixed infection of *P. falciparum* and other species.

Limitation and interferences

- The following anticoagulants have been validated for use with this test: heparin, EDTA and citrate.
- This test is designed to detect HRP2 and pLDH antigens of Malaria *P. falciparum* and *P. vivax*. A definitive clinical diagnosis should not be made based on the result of this test, but should only be made by a qualified physician after all clinical and laboratory findings have been evaluated.
- The prozone effect may cause false-negative result.
- Positive result with faint test line or false-negative is possible due to low parasite density.
- This test may produce a positive result after successful anti-malarial treatment. Therefore its use is not recommended for monitoring response to anti-malarial treatment.
- If specimens cannot be tested immediately, they should be refrigerated at 2-8°C for up to 3 days.



Rapydtest®

CareUS™ Malaria

Performance Characteristics

The careUS™ Malaria Combo Pf/PAN (HRP2/pLDH) Ag has been tested with positive and negative clinical samples conformed by microscopic examination.

Precision

Within- and between-run precisions have been evaluated by the testing 10 replicates of three specimens: a negative, a low positive and a strong positive. The agreement between the test results and the expected results were 100%.

Specimen	Positive	Negative	Sensitivity
P.falciparum	98	2	98%
P.vivax	96	4	96%

Specimen	Positive	Negative	Specificity
Negative	5	195	97.5%

95% CI, (97.05, 100) for P. falciparum positive result, (92.16, 99.84) for P. vivax positive results and (94.64, 99.36) for negative result.

References

- Valecha N., Eapen A., Devi C. Usha, Ravindran J., Aggarwal A., and Subbarao S. K. (2002). Field evaluation of the ICT Malaria P.f./P.v. immunochromatographic test in India. *Annals of Tropical Medicine & Parasitology*. 96: 333-336
- Iqbal J., Hira P. R., Sher A., and AL-Enezi A. A. 2001. Diagnosis of imported Malaria by Plasmodium Lactate Dehydrogenase (pLDH) and Histidine-Rich Protein 2 (PfHRP2)-based immunocapture assays. *American Journal of Tropical Medicine and Hygiene*. 64: 20-23
- Tjitra E., Suprianto S., Dyer M., Currie B. J. and Anstey N.M. (1999). Detection of Histidine-rich Protein 2 and panmalarial ICT MALARIA P.f./P.v. test antigens after chloroquine treatment of uncomplicated falciparum malaria does not reliably predict treatment outcome in eastern Indonesia. *Journal of Clinical microbiology*. 37: 2412-2417
- Panton L. J., PcPhie P., Maloy W. L., Welles T. E., Taylor D. W. and Howard R. J. (1989). Purification and partial characterization of an unusual protein of Plasmodium falciparum: histidine-rich protein II. *Molecular and Biochemical Parasitology*. 35: 149-160
- Leonard K. Basco, Frederique Marquet, Michael M.Makler, and Jacques Le Bras. : Plasmodium falciparum and Plasmodium vivax : Lactate Dehydrogenase Activity and its Application for in vitro Drug Susceptibility Assay. *Experimental Parasitology* 80, 260-271 (1995).

Ordering Information

PRODUCT	PACK SIZE	CODE
CareUS™ Malaria Rapydtest®, 30 individually wrapped strips, complete with Assay Buffer, Lancet, Pipette and Alcohol Swabs	30	1629
CareUS™ Malaria Rapydtest®, 5 individually wrapped strips, complete with Assay Buffer, Lancet, Pipette and Alcohol Swabs	5	16295

Products can be ordered direct from Apacor or from an appointed distributor


Visit our website for all the latest information www.apacor.com or e-mail on: orders@apacor.com

Description of Symbol Used

- Single use
- Expiry Date
- Batch code
- Manufacturer
- Catalogue number
- Keep dry
- CE mark
- Quantity in Box
- In-Vitro diagnostic medical device
- Storage Temperature limitation
- Consult instructions for use
- Do not use if package is damaged
- Keep away from sunlight
- Authorized representative in EU

WELLS BIO, INC
16, MAGOKJUNGANG 8-RO 1-GIL,
GANGSEO-GU, SEOUL, 07795
REPUBLIC OF KOREA

APACOR
 UNIT 5, SAPPHIRE CENTRE,
 FISHPONDS ROAD, WOKINGHAM,
 BERKSHIRE, RG41 2QL, UK
 TEL: +44 (0)118 979 5566
 FAX: +44 (0)118 979 5186



MDSS GmbH
Schiffaraben 41
30175 Hanover
Germany