



FOR FAECAL CONCENTRATION OF

HELMINTH OVA AND LARVAE / PROTOZOA CYSTS AND OOCYSTS

APACOR

Mini Parasep® SF
FAECAL PARASITE CONCENTRATOR



SELF STANDING
SAMPLE CHAMBER

MIXING CHAMBER

INTEGRAL
SPOON

Filter

A triple stage matrix filtration. Large particles are rejected without obscuring filtration. Recovery rate with Parasep® is comparable to traditional sieve method, ie: Ridley-Allen. The vertical filter enclosed design is patented.

Debris Trap

Rejected particles are trapped to prevent extrusion into the Sedimentation Cone during centrifugation.

Air/Liquid Seal and Safety Lock

The 'seal' prevents the release of biohazardous material. The 'lock' ensures the Mixing Chamber and Filter are removed together for safe disposal.

Sedimentation Cone

Sediment forms in the base of the cone allowing examination for the presence of helminth eggs or larvae and protozoa cysts or oocysts.

Health and Safety Benefits

- Totally enclosed/sealed process
- Reduced reagent volumes
- No cleaning required
- Single use, no sample contamination
- Ready to use systems available

Performance Benefits

- Optimum sample recovery
- Enhanced sample clarity
- Rapid four step process
- Human resources optimised
- Easy patient identification
- Fits all 15ml centrifuge buckets



PARASITOLOGY

SINGLE USE IN VITRO DIAGNOSTIC DEVICE



FAECAL PARASITE CONCENTRATOR
Mini Parasep® SF

Procedure

STEP 1 - SAMPLE PREPARATION

Preserved Samples

Shake or vortex the incoming preserved sample to thoroughly mix.

Transfer either 3ml of the emulsified stool into the Mini Parasep® SF mixing chamber.

In the event of:

Thick Stool Samples—please add 10 drops of Apacor Triton X solution, then please enclose and vortex/shake to emulsify prior to transferring the sample;

Liquid Stool Samples—please add 4ml instead of 3ml to ensure that a sediment is formed after centrifugation is performed.



STEP 2 - EMULSIFICATION

Seal the Mini Parasep® SF by screwing in the filter/sedimentation cone unit.



STEP 3 - CENTRIFUGATION


Invert Mini Parasep® SF and perform centrifugation at:

400g for two minutes.

NOTE: TO CALCULATE THE REQUIRED RPM FOR ANY CENTRIFUGE.

$$RPM = \sqrt{\frac{g}{1.12r}} \times 1000$$

RPM - rotor speed in revs/min.
 g - centrifugal force (max. 1000g)
 r - radius, horizontal distance between sedimentation cone tip and spindle centre measured in mm.




STEP 4 - EXAMINATION

Unscrew and discard the filter and mixing tube.

Decant the supernatant.

Transfer sediment to slide to perform examination where all temporary (Lugol's Iodine) and permanent staining techniques (Trichrome) can be conducted providing a universal fixative is used.

If a universal fixative is used in conjunction with this device, then molecular and EIA tests can be conducted from the sediment.



See label for storage conditions and expiry date. Please adhere to the following guidelines when handling Mini Parasep® SF. To avoid cross contamination the Mini Parasep® SF device should remain closed at all times except when introducing the sample or when retrieving the final concentrated sample for examination.

Mini Parasep® SF is available empty or reagent ready
Please ask for details

Products can be ordered direct from Apacor or from an appointed distributor
 Visit our website for all the latest information www.apacor.com or email on: orders@apacor.com

APACOR UNIT 5 SAPPHIRE CENTRE
 FISHPONDS ROAD, WOKINGHAM
 BERKSHIRE, RG41 2QL, UNITED KINGDOM
 TEL: +44 (0)118 979 5566
 FAX: +44 (0)118 979 5186

