


# Faecal Occult Blood (FOB) Rapydtest®

For the detection of human  
haemoglobin in faeces

APACOR



MICROBIOLOGY  
SINGLE USE IN VITRO  
DIAGNOSTIC DEVICE

# Faecal Occult Blood (FOB)

Rapydtest®

## Intended use

The FOB Rapydtest® is an immunochemical device for the qualitative detection of faecal occult blood. It is a useful aid to detect bleeding caused by a number of gastrointestinal disorders and is recommended for use in routine physical examinations, routine hospital testing and screening for colorectal cancer or gastrointestinal bleeding from any source.






## Performance characteristics

A total of 251 samples were collected from two hospitals and tested by the FOB Rapydtest® and by a leading commercial fob rapid test.

**Relative sensitivity - 100%**

**Relative specificity - 99.4%**

**Overall agreement - 99.6%**

	<b>Negative result</b>
	<b>Positive result</b>
	
	<b>Invalid result</b>
	





\*\*For liquid samples



2 drops  
70-90ul



2 drops



2 drops



10 minutes

Code	Description	Pack size
1642	Faecal Occult Blood (FOB) Rapydtest®	25

Visit [www.apacor.com](http://www.apacor.com) for all our latest information or contact us at [sales@apacor.com](mailto:sales@apacor.com)

[www.apacor.com](http://www.apacor.com)

APA-BR-020/04



SAPPHIRE CENTRE  
FISHPONDS ROAD, WOKINGHAM  
BERKSHIRE, RG41 2QL, ENGLAND

TEL: +44 (0)118 979 5566

FAX: +44 (0)118 979 5186

APACOR