



Faecal Occult Blood (FOB) Rapydtest®

For the detection of human
haemoglobin in faeces

APACOR



MICROBIOLOGY
SINGLE USE IN VITRO
DIAGNOSTIC DEVICE

Faecal Occult Blood (FOB)

Rapydtest®






Intended use

The FOB Rapydtest® is an immunochemical device for the qualitative detection of faecal occult blood. It is a useful aid to detect bleeding caused by a number of gastrointestinal disorders and is recommended for use in routine physical examinations, routine hospital testing and screening for colorectal cancer or gastrointestinal bleeding from any source.

Performance characteristics

A total of 251 samples were collected from two hospitals and tested by the FOB Rapydtest® and by a leading commercial fob rapid test.

Relative sensitivity - 100%
Relative specificity - 99.4%
Overall agreement - 99.6%

	Negative result
	Positive result
	
	Invalid result
	




**For liquid samples



2 drops
70-90ul



2 drops



2 drops



10 minutes

Code	Description	Pack size
1642	Faecal Occult Blood (FOB) Rapydtest®	25

Visit www.apacor.com for all our latest information or contact us at sales@apacor.com