

30ml Transportation Vials

for Transportation of faecal specimens
for intestinal parasite examination

APACOR



PARASITOLOGY

30ml Transportation Vials





for Transportation of faecal specimens

In the laboratory, intestinal parasites are confirmed by the identification of cysts, protozoan trophozoites, larvae and helminth eggs.

The priorities of these clinical laboratories do not always permit the immediate examination of a fresh sample. The prompt collection and Transportation of these samples cannot always be guaranteed.

Furthermore, freezing, refrigeration and/or incubation of samples cannot ensure the full recovery of parasites at all stages of identification.

The Apacor Transportation vial will preserve the intestinal parasites in the faecal material until such time it can be examined by a qualified parasitologist.



Visit www.apacor.com for all our latest information or contact us at sales@apacor.com

Code	Description	Pack size
149970	30ml Transportation vial (Clean)	40
148998	30ml Transportation vial with 15ml Formalin & Triton X	40
249400	30ml Transportation vial with 15ml SAF	40
249420	30ml Transportation vial with 15ml AlcorFix™	40
160001	30ml Transportation vial with 15ml Apafix™	40